

Anti-Birc7 Polyclonal Antibody For Wb, Elisa - Q96Ca5

Item Number: CM0026285



Introduction

BIRC7 Rabbit pAb validated for WB and ELISA detection of human and mouse BIRC7 (Livin/ML-IAP). Ideal for apoptosis and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	BIRC7 Rabbit pAb
Remarks/Alias	KIAP; LIVIN; MLIAP; RNFS0; ML-IAP; BIRC7
Species	Human
Gene ID (Human)	79444
Gene ID	79444
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human BIRC7 (NP_647478.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96CA5
Protein Weight	33kDa
Shipping	Ice bag