


Anti-Blbp/Fabp7 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - O15540

Item Number: CM0006688



Introduction

High-specificity rabbit monoclonal antibody targeting human BLBP/FABP7, a brain lipid-binding protein involved in CNS development. Validated for WB, IF-P, and ELISA; cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	BLBP/FABP7 Rabbit mAb
Remarks/Alias	MRG; BLBP; FABPB; B-FABP; BLBP/FABP7
Species	Human
GeneID (Human)	2173
GeneID	2173
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human BLBP/FABP7 (NP_001437.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15540
Protein Weight	14KDa/18KDa
Shipping	Ice bag