

Anti-Bloc1S2 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Qny1

Item Number: CM0028814

Introduction

Rabbit polyclonal antibody against human BLOC1S2, validated for WB and ELISA in human, mouse, and rat. Ideal for studying lysosome-related organelle biogenesis, neurite extension, and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	BLOC1S2 Rabbit pAb
Remarks/Alias	CEAP, BLOS2, BORCS2, CEAP11, BLOC1S2
Species	Human
GeneID (Human)	282991
GeneID	282991
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 44-142 of human BLOC1S2 (NP_776170.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6QNY1
Protein Weight	16kDa
Shipping	Ice bag