

Anti-Bloc1S3 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q6Qny0

Item Number: CM0023339

Introduction

Rabbit polyclonal antibody against human BLOC1S3, validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for lysosome-related organelle biogenesis research, including BLOC-1 complex studies.

[Learn More](#)

Parameter	Value
Product Name	BLOC1S3 Rabbit pAb
Remarks/Alias	RP; HPS8; BLOS3; BLOC1S3
Species	Human
GeneID (Human)	388552
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-202 of human BLOC1S3 (NP_997715.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6QNY0
Protein Weight	21kDa
Shipping	Ice bag