

Anti-Bloc1S6 Rabbit Polyclonal Antibody For Wb And Elisa - Q9UI45

Item Number: CM0026193

Introduction

BLOC1S6 Rabbit Polyclonal Antibody validated for Western Blot (WB) and ELISA. Detects human, mouse, and rat BLOC1S6/pallidin, a subunit of the BLOC-1 complex essential for lysosome-related organelle biogenesis and intracellular trafficking.

[Learn More](#)

Parameter	Value
Product Name	BLOC1S6 Rabbit pAb
Remarks/Alias	PA; HPS9; PLDN; BLOS6; PALLID; BLOC1S6
Species	Human
GeneID	26258
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human BLOC1S6 (NP_036520.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UL45
Protein Weight	20kDa
Shipping	Ice bag