


Anti-Blvrb Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P30043

Item Number: CM0025745



Introduction

Rabbit polyclonal antibody against human BLVRB (FLR, BVRB), validated in WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying heme degradation, redox biology, and insulin signaling.

[Learn More](#)

Parameter	Value
Product Name	BLVRB Rabbit pAb
Remarks/Alias	FLR; BVRB; SDR43U1; HEL-S-10; BLVRB
Species	Human
Gene ID (Human)	645
Gene ID	645
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human BLVRB (NP_000704.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30043
Protein Weight	22 kDa
Shipping	Ice bag