

Anti-Bmal1 Monoclonal Antibody For Wb, Elisa - 000327

Item Number: CM0005507



Introduction

High-quality rabbit monoclonal antibody targeting human BMAL1 (gene ARNTL), validated for WB and ELISA. Detects endogenous BMAL1 across human, mouse, and rat species. Key reagent for circadian rhythm, metabolism, and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	BMAL1 Rabbit mAb
Remarks/Alias	TIC; JAP3; MOP3; ARNTL; PASD3; ARNTL1; BMAL1c; bHLHe5; BMAL1
Species	Human
GeneID (Human)	406
GeneID	406
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 500-600 of human BMAL1 (000327))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	000327
Protein Weight	69kDa
Shipping	Ice bag