

Anti-Bmal1 Polyclonal Antibody For Wb, Elisa - 000327

Item Number: CM0022082



Introduction

High-quality anti-BMAL1 polyclonal antibody (Rabbit) validated for WB and ELISA, reacting with human, mouse, rat. Ideal for circadian rhythm research and IVD applications.

[Learn More](#)

Parameter	Value
Product Name	BMAL1 Rabbit pAb
Remarks/Alias	TIC; JAP3; MOP3; ARNTL; PASD3; ARNTL1; BMAL1c; bHLHe5; BMAL1
Species	Human
GeneID (Human)	406
GeneID	406
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-510 of human BMAL1 (NP_001284648.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00327
Protein Weight	69kDa
Shipping	Ice bag