


Anti-Bmi1 Polyclonal Antibody For Wb, If/Icc, Elisa - P35226

Item Number: CM0015006



Introduction

Rabbit polyclonal antibody against human Bmi1 (PCGF4/RNF51), suitable for WB, IF/ICC, ELISA. Target polycomb complex protein BMI-1, a chromatin regulator involved in gene repression and stem cell maintenance.

[Learn More](#)

Parameter	Value
Product Name	Bmi1 Rabbit pAb
Remarks/Alias	PCGF4; RNF51; FLV12/BMI1; flvi-2/bmi-1; Bmi1
Species	Human
Gene ID (Human)	648
Gene ID	648
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 120-326 of human Bmi1 (NP_005171.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P35226
Protein Weight	37kDa
Shipping	Ice bag