


Anti-Bmi1 Polyclonal Antibody For Wb, Ip, Elisa - P35226

Item Number: CM0015782



Introduction

Rabbit polyclonal antibody targeting human Bmi1 (PCGF4/RNF51), a core Polycomb group protein. Validated for WB, IP, ELISA; cross-reacts with human, mouse, rat. Ideal for epigenetic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Bmi1 Rabbit pAb
Remarks/Alias	PCGF4; RNF51; FLV12/BMI1; flvi-2/bmi-1; Bmi1
Species	Human
Gene ID	648
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Bmi1 (NP_005171.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35226
Protein Weight	37kDa
Shipping	Ice bag