


Anti-Bmi1 Rabbit Monoclonal Antibody For Wb, Elisa - P35226

Item Number: CM0002611



Introduction

Rabbit monoclonal antibody against human Polycomb complex protein BMI-1 (Bmi1, PCGF4). Validated for WB and ELISA, reacts with human and rat. ~37kDa, key for chromatin remodeling and gene silencing research.

[Learn More](#)

Parameter	Value
Product Name	Bmi1 Rabbit mAb
Remarks/Alias	PCGF4, RNF51, FLVI2/BMI1, flvi-2/bmi-1, Bmi1
Species	Human
GeneID (Human)	648
GeneID	648
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 227-326 of human Bmi1 (P35226)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P35226
Protein Weight	37kDa
Shipping	Ice bag