


Anti-Bmp1 Polyclonal Antibody For If/Icc, Elisa - P13497

Item Number: CM0019263



Introduction

BMP1 rabbit polyclonal antibody (pAb) for human samples, validated in IF/ICC and ELISA. Targets procollagen C-proteinase involved in ECM formation and bone repair. Suitable for matrix biology and developmental research.

[Learn More](#)

Parameter	Specification
Product Name	BMP1 Rabbit pAb
Remarks/Alias	PCP; TLD; OI13; PCP2; PCOLC; BMP1
Species	Human
GeneID (Human)	649
GeneID	649
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 725-986 of human BMP1 (NP_006120.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P13497
Protein Weight	111kDa
Shipping	Ice bag