

Anti-Bmp10 Polyclonal Antibody For Wb, Elisa - 095393

Item Number: CM0017009



Introduction

Rabbit polyclonal antibody targeting human BMP10, validated for Western blot and ELISA. Cross-reacts with rat. Ideal for investigating cardiac development, endothelial cell migration, and BMP signaling. Recombinant immunogen, predicted 48kDa.

[Learn More](#)

Parameter	Value
Product Name	BMP10 Rabbit pAb
Remarks/Alias	BMP10
Species	Human
GeneID (Human)	27302
GeneID	27302
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-316 of human BMP10 (NP_055297.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	O95393
Protein Weight	48kDa
Shipping	Ice bag