

Anti-Bmp15 Rabbit Monoclonal Antibody For Wb, Elisa - O95972

Item Number: CM0010729

Introduction

Recombinant rabbit monoclonal antibody against human BMP15 (O95972). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for folliculogenesis and granulosa cell growth research.

[Learn More](#)

Parameter	Value
Product Name	BMP15 Rabbit mAb
Remarks/Alias	ODG2; POF4; GDF9B; BMP15
Species	Human
GeneID (Human)	9210
GeneID	9210
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 293-392 of human BMP15 (O95972).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95972
Protein Weight	45kDa
Shipping	Ice bag