

Anti-Bmp2 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P12643

Item Number: CM0010144



Introduction

Rabbit monoclonal antibody targeting human BMP2 for WB, IF/ICC, IP, ELISA. Cross-reactivity with mouse and rat. Recombinant protein immunogen. Suitable for bone morphogenesis, osteogenesis, and TGF-beta signaling studies.

[Learn More](#)

Parameter	Value
Product Name	BMP2 Rabbit mAb
Remarks/Alias	BDA2; BMP2A; SSFSC; SSFSC1; BMP2
Species	Human
GeneID (Human)	650
GeneID	650
Immunogen	Recombinant protein Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12643
Protein Weight	45kDa
Shipping	Ice bag