

Anti-Bmp2 Rabbit Polyclonal Antibody For Wb, Elisa - P12643

Item Number: CM0018057



Introduction

Rabbit polyclonal antibody against human BMP2, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. A key research tool for investigating osteogenesis, bone morphogenetic protein signaling, and developmental processes.

[Learn More](#)

Parameter	Description
Product Name	BMP2 Rabbit pAb
Remarks/Alias	BDA2; BMP2A; SSFSC; SSFSC1; BMP2
Species	Human
GeneID (Human)	650
GeneID	650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human BMP2 (NP_001191.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12643
Protein Weight	45kDa
Shipping	Ice bag