

Anti-Bmp4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P12644

Item Number: CM0013390

Introduction

Rabbit polyclonal antibody targeting human BMP4, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying TGF-beta signaling, developmental biology, and osteogenesis.

[Learn More](#)

Parameter	Value
Product Name	BMP4 Rabbit pAb
Remarks/Alias	ZYME; BMP2B; OFC11; BMP2B1; MCOP56; BMP4
Species	Human
GeneID (Human)	652
GeneID	652
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human BMP4 (NP_001193.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12644
Protein Weight	47kDa
Shipping	Ice bag