

Anti-Bmp4 Rabbit Monoclonal Antibody For Wb And Elisa - P12644

Item Number: CM0002434

Introduction

BMP4 rabbit monoclonal antibody validated for WB and ELISA in Human, Mouse, Rat. Bone morphogenetic protein 4 is a TGF-beta superfamily member involved in development and osteogenesis. Available from CamelBio as a one-stop IVD raw material partner.

[Learn More](#)

Parameter	Value
Product Name	BMP4 Rabbit mAb
Remarks/Alias	ZYME; BMP2B; OFC11; BMP2B1; MCOPS6; BMP4
Species	Human
Gene ID (Human)	652
Gene ID	652
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human BMP4 (NP_001193.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12644
Protein Weight	47kDa
Shipping	Ice bag