

Anti-Bmp4 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P12644

Item Number: CM0003433

Introduction

BMP4 rabbit monoclonal antibody validated for Western blot, IHC (paraffin), IF/ICC, and ELISA; specifically detects human BMP4, a TGF-beta growth factor essential for development, osteogenesis, and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	BMP4 Rabbit mAb
Remarks/Alias	ZYME, BMP2B, OFC11, BMP2B1, MCOPS6, BMP4
Species	Human
GeneID (Human)	652
GeneID	652
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 309-408 of human BMP4 (P12644).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P12644
Protein Weight	47kDa
Shipping	Ice bag