

Anti-Bmp5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22003

Item Number: CM0029203

Introduction

High-quality rabbit polyclonal antibody against human BMP5, validated for WB, IHC-P, and ELISA. Detects BMP5 in human, mouse, and rat samples. Ideal for bone morphogenetic protein signaling research.

[Learn More](#)

Parameter	Value
Product Name	BMP5 Rabbit pAb
Remarks/Alias	BMP5
Species	Human
Gene ID (Human)	653
Gene ID	653
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 317-454 of human BMP5 (NP_066551.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22003
Protein Weight	52kDa
Shipping	Ice bag