


Anti-Bmp6 Monoclonal Antibody For Wb, If/Icc, Elisa - P22004

Item Number: CM0012318



Introduction

High-quality anti-BMP6 rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for bone development, iron metabolism, and TGF-beta signaling research.

[Learn More](#)

Parameter	Value
Product Name	BMP6 Rabbit mAb
Remarks/Alias	IO; VGR; VGR1; BMP6
Species	Human
GeneID (Human)	654
GeneID	654
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22004
Protein Weight	57kDa
Shipping	Ice bag