


Anti-Bmp6 Polyclonal Antibody For If/Icc And Elisa - P22004

Item Number: CM0018568



Introduction

Rabbit polyclonal antibody against human BMP6, validated for immunofluorescence (IF/ICC) and ELISA. Reacts with human and mouse BMP6 (57 kDa). Ideal for bone development and iron metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	BMP6 Rabbit pAb
Remarks/Alias	IO; VGR; VGR1; BMP6
Species	Human
GeneID (Human)	654
GeneID	654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 375-513 of human BMP6 (NP_001709.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P22004
Protein Weight	57kDa
Shipping	Ice bag