


Anti-Bmp6 Polyclonal Antibody For Wb And Elisa - P22004

Item Number: CM0025938



Introduction

High-quality rabbit polyclonal antibody against human BMP6, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for bone development, iron metabolism, and TGF-beta signaling research.

[Learn More](#)

Parameter	Value
Product Name	BMP6 Rabbit pAb
Remarks/Alias	IO; VGR; VGR1; BMP6
Species	Human
GeneID (Human)	654
GeneID	654
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 375-513 of human BMP6 (NP_001709.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22004
Protein Weight	57kDa
Shipping	Ice bag