


Anti-Bmp7 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P18075

Item Number: CM0017104



Introduction

Rabbit polyclonal antibody against BMP7, validated for WB, IHC-P, ELISA in human, mouse, rat samples. Suitable for studying TGF-beta signaling, development, and osteogenesis.

[Learn More](#)

Parameter	Value
Product Name	BMP7 Rabbit pAb
Remarks/Alias	OP-1; BMP7
Species	Human
GeneID (Human)	655
GeneID	655
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-251 of BMP7 (NP_001710.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18075
Protein Weight	49kDa
Shipping	Ice bag