

Anti-Bmpr1A Rabbit Polyclonal Antibody For Wb, Elisa - P36894

Item Number: CM0022602



Introduction

BMPR1A (ALK3/CD292) rabbit polyclonal antibody for WB and ELISA applications, reactive in human, mouse, and rat. Useful for investigating BMP receptor signaling, skeletal development, and bone morphogenetic protein pathways.

[Learn More](#)

Parameter	Value
Product Name	BMPR1A Rabbit pAb
Remarks/Alias	ALK3; SKR5; ALK-3; CD292; ACVRLK3; BMPR-1A; 10q23del; BMPR1A
Species	Human
GeneID (Human)	657
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-200 of human BMPR1A (NP_004320.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36894
Protein Weight	60kDa
Shipping	Ice bag