

Anti-Bmpr2 Rabbit Polyclonal Antibody - Q13873

Item Number: CM0021602



Introduction

High-quality rabbit polyclonal anti-BMPR2 antibody for research use. Validated in WB, IF-P, IHC-P, ELISA. Recognizes human, mouse, rat BMP type II receptor (BMPR-II). Recombinant immunogen (aa27-150).

[Learn More](#)

Parameter	Value
Product Name	BMPR2 Rabbit pAb
Remarks/Alias	BMR2; PPH1; BMPR3; BRK-3; POVD1; T-ALK; BMPR-II; BMPR2
Species	Human
GeneID (Human)	659
GeneID	659
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-150 of human BMPR2 (NP_001195.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13873
Protein Weight	115kDa
Shipping	Ice bag