


Anti-Bmx Polyclonal Antibody For Wb, Elisa - P51813

Item Number: CM0025857



Introduction

Rabbit polyclonal antibody against BMX (ETK, PSCTK2, PSCTK3), validated for WB and ELISA, cross-reacts with human and mouse. Useful for studying signal transduction, angiogenesis, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	BMX Rabbit pAb
Remarks/Alias	ETK; PSCTK2; PSCTK3; BMX
Species	Human
GeneID (human)	660
GeneID	660
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human BMX (NP_001712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51813
Protein Weight	78kDa
Shipping	Ice bag