

Anti-Bnip1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q12981

Item Number: CM0012748

Introduction

Rabbit monoclonal antibody targeting human BNIP1 (SEC20) for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for research on ER stress, apoptosis, and mitochondrial autophagy.

[Learn More](#)

Parameter	Details
Product Name	BNIP1 Rabbit mAb
Remarks/Alias	NIP1; SEC20; TRG-8; BNIP1
Species	Human
GeneID (Human)	662
GeneID	662
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human BNIP1 (Q12981)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q12981
Protein Weight	26kDa
Shipping	Ice bag