

Anti-Bnip1 Rabbit Polyclonal Antibody For Wb And Elisa - Q12981

Item Number: CM0028522

Introduction

Rabbit polyclonal antibody targeting human BNIP1 (SEC20/NIP1) protein. Suitable for Western blot (WB) and ELISA applications, with cross-reactivity against mouse and rat. Supports research into endoplasmic reticulum membrane fusion, apoptosis, and mitochondrial autophagy.

[Learn More](#)

Parameter	Value
Product Name	BNIP1 Rabbit pAb
Remarks/Alias	NIP1; SEC20; TRG-8; BNIP1
Species	Human
GeneID (human)	662
GeneID	662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human BNIP1 (NP_001196.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12981
Protein Weight	26kDa
Shipping	Ice bag