

Anti-Bnip2 Polyclonal Antibody For Wb & Elisa - Q12982

Item Number: CM0024784

Introduction

Rabbit polyclonal antibody targeting human BNIP2, a BCL-2/E1B 19 kDa-interacting protein involved in apoptosis suppression. Validated for Western blot (WB) and ELISA, with cross-reactivity to human and mouse. Predicted molecular weight 36 kDa.

[Learn More](#)

Parameter	Value
Product Name	BNIP2 Rabbit pAb
Remarks/Alias	NIP2; BNIP-2; BNIP2
Species	Human
Gene ID (Human)	663
Gene ID	663
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 125-280 of human BNIP2 (NP_004321.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12982
Protein Weight	36kDa
Shipping	Ice bag