


Anti-Bnip3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12983

Item Number: CM0030765



Introduction

Rabbit polyclonal antibody targeting human BNIP3, validated for WB, IHC-P, and ELISA applications. Reacts with human, mouse, and rat. Suitable for apoptosis, autophagy, and mitochondrial quality control studies.

[Learn More](#)

Parameter	Value
Product Name	BNIP3 Rabbit pAb
Remarks/Alias	NIP3; HABON; BNIP3
Species	Human
GeneID (human)	664
GeneID	664
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-120 of human BNIP3 (NP_004043.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12983
Protein Weight	22kDa
Shipping	Ice bag