

Anti-Bnip3 Rabbit Polyclonal Antibody For If/Icc And Elisa - Q12983

Item Number: CM0015112



Introduction

High-quality rabbit polyclonal antibody targeting human BNIP3, validated for IF/ICC and ELISA. Ideal for investigating apoptosis, mitochondrial quality control, and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	BNIP3 Rabbit pAb
Remarks/Alias	NIP3; HABON; BNIP3
Species	Human
GeneID (Human)	664
GeneID	664
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human BNIP3 (NP_004043.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q12983
Protein Weight	22kDa
Shipping	Ice bag