

## Anti-Bnip3L Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60238

Item Number: CM0028346

### Introduction

Rabbit polyclonal antibody against human BNIP3L (NIX). Validated for WB, IHC-P, IF/ICC, ELISA. Detects BNIP3L across human, mouse, and rat samples. Ideal for apoptosis and mitochondrial quality control studies.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | BNIP3L Rabbit pAb                                                                                                  |
| Remarks/Alias    | NIX; BNIP3a; BNIP3L                                                                                                |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 665                                                                                                                |
| GeneID           | 665                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 95-180 of human BNIP3L (NP_004322.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | O60238                                                                                                             |
| Protein Weight   | 24kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |