

Anti-Bnp Polyclonal Antibody For Wb, If, Ihc, Elisa - P16860

Item Number: CM0020551



Introduction

BNP Rabbit Polyclonal Antibody detects the cardiac hormone BNP in human, mouse, and rat samples. Validated for WB, IF/ICC, IF-P, IHC-P, and ELISA. Ideal for cardiovascular and renal research.

[Learn More](#)

Parameter	Value
Product Name	BNP Rabbit pAb
Remarks/Alias	BNP; Iso-ANP
Species	Human
GeneID (Human)	4879
GeneID	4879
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16860
Protein Weight	15 kDa
Shipping	Ice bag