

Anti-Bnp Rabbit Monoclonal Antibody For Wb, Elisa - P16860

Item Number: CM0007460



Introduction

Rabbit monoclonal antibody targeting human BNP (Iso-ANP), validated for WB and ELISA applications. Cross-reacts with mouse and rat. Ideal for cardio-renal and natriuretic peptide research.

[Learn More](#)

Parameter	Value
Product Name	BNP Rabbit mAb
Remarks/Alias	BNP; Iso-ANP
Species	Human
Gene ID (Human)	4879
Gene ID	4879
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 35-134 of human BNP (NP_002512.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16860
Protein Weight	15kDa
Shipping	Ice bag