

Anti-Bob1 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q16633

Item Number: CM0010519



Introduction

Rabbit monoclonal antibody targeting human BOB1 (POU2AF1), validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for B-cell transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	BOB1 Rabbit mAb
Remarks/Alias	BOB1; OBF1; OCAB; OBF-1
Species	Human
Gene ID (Human)	5450
Gene ID	5450
Immunogen	Synthetic peptide corresponding to amino acids 150-250 of human BOB1 (NP_006226.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16633
Protein Weight	27kDa
Shipping	Ice bag