


Anti-Bok Rabbit Monoclonal Antibody For Wb And Elisa - Q9Umx3

Item Number: CM0006990



Introduction

High-specificity rabbit monoclonal antibody targeting human BOK (Bcl-2-related ovarian killer protein). Validated for WB and ELISA with cross-reactivity across human, mouse, and rat. Ideal for apoptosis and mitochondrial pathway studies.

[Learn More](#)

Parameter	Value
Product Name	BOK Rabbit mAb
Remarks/Alias	BOKL; BCL2L9; BOK
Species	Human
GeneID (human)	666
GeneID	666
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UMX3
Protein Weight	23kDa
Shipping	Ice bag