


Anti-Bok Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Umx3

Item Number: CM0020643



Introduction

Rabbit polyclonal antibody targeting human BOK (BCL2L9). Applicable in WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Supports apoptosis and oncology research.

[Learn More](#)

Parameter	Value
Product Name	BOK Rabbit pAb
Remarks/Alias	BOKL; BCL2L9; BOK
Species	Human
Gene ID (Human)	666
Gene ID	666
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human BOK (NP_115904.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UMX3
Protein Weight	23 kDa
Shipping	Ice bag