


Anti-Bola3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q53S33

Item Number: CM0024332



Introduction

Rabbit polyclonal antibody against human BOLA3, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for mitochondrial Fe-S cluster assembly studies and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	BOLA3 Rabbit pAb
Remarks/Alias	MMD52; BOLA3
Species	Human
GeneID (Human)	388962
GeneID	388962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human BOLA3 (NP_997717.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53S33
Protein Weight	12kDa
Shipping	Ice bag