

Anti-Boll Polyclonal Antibody For Wb, Elisa - Q8N9W6

Item Number: CM0027803

Introduction

Rabbit polyclonal antibody against human BOLL (BOULE), validated for Western blot and ELISA. Detects BOLL in human, mouse, and rat samples. Suitable for spermatogenesis and RNA-binding protein research.

[Learn More](#)

Parameter	Specification
Product Name	BOLL Rabbit pAb
Remarks/Alias	BOULE; BOLL
Species	Human
Gene ID (Human)	66037
Gene ID	66037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human BOLL (NP_149019.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N9W6
Protein Weight	31kDa
Shipping	Ice bag