

Anti-Boule/Boll Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q8N9W6

Item Number: CM0012261

Introduction

Recombinant rabbit monoclonal antibody against human BOLL (BOULE), validated for WB, IP, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for spermatogenesis and fertility research.

[Learn More](#)

Parameter	Value
Product Name	BOULE/BOLL Rabbit mAb
Remarks/Alias	BOULE
Species	Human
GeneID (human)	66037
GeneID	66037
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human BOLL (NP_932074.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N9W6
Protein Weight	31kDa
Shipping	Ice bag