


Anti-Bpgm Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P07738

Item Number: CM0029751



Introduction

Rabbit polyclonal antibody targeting human BPGM (bisphosphoglycerate mutase). Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for studying oxygen affinity regulation and erythrocyte disorders.

[Learn More](#)

Parameter	Specification
Product Name	BPGM Rabbit pAb
Remarks/Alias	DPGM; ECYT8; BPGM
Species	Human
GeneID (human)	669
GeneID	669
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-259 of human BPGM (NP_001715.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07738
Protein Weight	30kDa
Shipping	Ice bag