


Anti-Bpi Rabbit Polyclonal Antibody For Wb, Elisa - P17213

Item Number: CM0029472



Introduction

Rabbit polyclonal antibody targeting human bactericidal permeability-increasing protein (BPI). Validated for WB and ELISA. Cross-reacts with mouse. Ideal for innate immunity and Gram-negative bacteria research.

[Learn More](#)

Parameter	Value
Product Name	BPI Rabbit pAb
Remarks/Alias	rBPI; BPIFD1; BPI
Species	Human
Gene ID	671
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 248-487 of human BPI (NP_001716.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17213
Protein Weight	54kDa
Shipping	Ice bag