

Anti-Bpifa1 Polyclonal Antibody For Wb, Elisa - Q9Np55

Item Number: CM0029009



Introduction

Rabbit polyclonal antibody targeting human BPIFA1 (SPLUNC1), validated for WB and ELISA. Cross-reacts with Mouse, Rat. Suitable for investigating innate immunity, airway biology, and antimicrobial roles.

[Learn More](#)

Parameter	Value
Product Name	BPIFA1 Rabbit pAb
Remarks/Alias	LUNX, NASG, PLUNC, SPURT, SPLUNC1, bA49G10.5, BPIFA1
Species	Human
Gene ID	51297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-256 of human BPIFA1 (NP_057667.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NP55
Protein Weight	27kDa
Shipping	Ice bag