

Anti-Bpifb2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N4F0

Item Number: CM0013244

Introduction

Rabbit polyclonal antibody targeting human BPIFB2 (RYSR/BPIL1/LPLUNC2) for Western blot and ELISA. Cross-reacts with mouse and rat. A secreted glycoprotein highly expressed in tonsils, suitable for innate immunity and epithelial biology research.

[Learn More](#)

Parameter	Value
Product Name	BPIFB2 Rabbit pAb
Remarks/Alias	RYSR; BPIL1; LPLUNC2; C20orf184; dj726C3.2; BPIFB2
Species	Human
Gene ID (Human)	80341
Gene ID	80341
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-340 of human BPIFB2 (NP_079503.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8N4F0
Protein Weight	49kDa
Shipping	Ice bag