

Anti-Brachyury Polyclonal Antibody For Wb, If/Icc, Elisa - O15178

Item Number: CM0022783



Introduction

High-quality Brachyury rabbit polyclonal antibody from CamelBio, validated for WB, IF/ICC, ELISA in human and rat. Essential tool for studying mesoderm formation and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	Brachyury Rabbit pAb
Remarks/Alias	T; TFT; SAVA; Brachyury
Species	Human
GeneID	6862
Immunogen	Recombinant Brachyury fusion protein containing a sequence corresponding to amino acids 226-435 of human Brachyury (NP_003172.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O15178
Protein Weight	47kDa
Shipping	Ice bag