

Anti-Brachyury Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O15178

Item Number: CM0008102



Introduction

Recombinant rabbit monoclonal antibody targeting human Brachyury (T-box transcription factor T), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for mesoderm differentiation and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Brachyury Rabbit mAb
Remarks/Alias	T; TFT; SAVA; Brachyury
Species	Human
Gene ID (Human)	6862
Gene ID	6862
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 266-435 of human Brachyury (O15178)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15178
Protein Weight	47kDa
Shipping	Ice bag