

Anti-Braf V600E Monoclonal Antibody For Wb - P15056

Item Number: CM0007013



Introduction

High-specificity rabbit monoclonal antibody targeting BRAF V600E (84 kDa). Suitable for WB, IF/ICC, ELISA; human reactive. Ideal for oncology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	BRAF V600E Rabbit mAb
Remarks/Alias	NS7; B-raf; BRAF1; RAFB1; B-RAF1; BRAF-1; BRAF V600E
Species	Human
GeneID (Human)	673
GeneID	673
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P15056
Protein Weight	84kDa
Shipping	Ice bag