

Anti-Brap (Imp/Brap2/Rnf52) Polyclonal Antibody For Wb, Ihc-P, Elisa - Q7Z569

Item Number: CM0023808

Introduction

Rabbit polyclonal antibody against human BRAP (BRCA1-associated protein, IMP, BRAP2, RNF52). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for Ras signaling, ubiquitination, and nuclear transport research.

[Learn More](#)

Parameter	Value
Product Name	BRAP Rabbit pAb
Remarks/Alias	IMP; BRAP2; RNF52; BRAP
Species	Human
Gene ID (Human)	8315
Gene ID	8315
Immunogen	Recombinant Protein - Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human BRAP (NP_006759.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z569
Protein Weight	67kDa
Shipping	Ice bag