

Anti-Brca1 Polyclonal Antibody For Wb And Elisa - P38398

Item Number: CM0014185



Introduction

Rabbit polyclonal antibody against human BRCA1, validated for WB and ELISA. Reacts with human and mouse. Immunogen: synthetic peptide within aa 100-200. For studying DNA repair and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	BRCA1 Rabbit pAb
Remarks/Alias	IRIS; PSCP; BRCAI; BRCC1; FANCS; PNCA4; RNF53; BROVCA1; PPP1R53; BRCA1
Species	Human
GeneID (Human)	672
GeneID	672
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human BRCA1 (NP_009225.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P38398
Protein Weight	208kDa
Shipping	Ice bag